

PROTOCOL FOR PI STAINING OF DNA

- 1) Fix cells in 4% paraformaldehyde with 0.05% glutaraldehyde in PBS for 15 minutes.
- 2) Wash in PBS with 0.1% triton for 5 minutes.
- 3) Wash in PBS 2 x 5 minutes.
- 4) RNA digestion: 400ug/mg RNase A at 37°C for 30 minutes.
- 5) Wash in PBS 2 x 5 minutes.
- 6) Stain with PI (0.25mg/ml for 30 minutes) at RT.
- 7) Wash with PBS 3 x 5 minutes.
- 8) Fix in 2% glutaraldehyde for 5-15 minutes.
- 9) Wash in cacodylate then photooxidize with DAB 1mg/ml.
- 10) Wash in cacodylate 3 x 5 minutes.
- 11) Postfix in 1% OsO₄ in cacodylate for 30 minutes.
- 12) Standard EM procedure.